

punctuality is an indication of competence


technical data of cryogenic tunnel freezer*

Model	Production Capacity	Width, Overall Belt "A"	Width, Usable Belt* "B"	Length, Usable Belt** "C"	Product Clearance "D"	Length, Overall*** "E"	Width, Overall "F"	Height, Overall Open "G"	Height, Inlet*** "H"	Height, Outlet*** "J"	Weight	Cooling Medium	Power Supply
	lbs/hr	in	in	ft	in	in	in	in	in	in			
LT24x20	4000	26	24	20	6.5	26	66	100	45	45	4000	LIN or LCO ₂	480V
LT36x20	6000	38	36	20	6.5	26	78	100	45	45	7000	LIN or LCO ₂	480V
LT48x20	8000	50	48	20	6.5	26	90	100	45	45	10000	LIN or LCO ₂	480V
LT58x20	10000	60	58	20	6.5	26	108	100	45	45	13000	LIN or LCO ₂	480V

* Values listed for single belt tunnels. Dual-belt options available

** Values listed for 20 ft length body design. Body sections come in 10' increments. Add approx. 6' for end blocks

*** Overall height may be altered slightly by adjusting legs



* US specifications and build are different from the European one